



Hard Gelatin Capsule Machine

7D2A 2-Deck Elongated Model



Hard Gelatin Capsule Machine

7D2A 2-Deck Elongated Model

Efficient Hard Gel Capsule Drying

The 7D2A 2-Deck TES elongated model has been designed to provide automated capsule filling capabilities for Hard Gel Capsules. Designed to provide high efficiency for production scale use meeting GMP standards for both pharmaceutical and nutraceutical industries.

The system consists of:

- Two deck machine: Upper Deck: 30 tables
 Lower Deck: 25 tables
- Control of all production and process parameters using high-tech control panel.
- Very easy and flexible operation using a Programmable Logic Controller (PLC) with touchscreen.

Features:

- Two sets of automatics.
- Dipping 7-bars.
- 28 bars on table.
- Pin bar requirements: 3,208 pieces for 7-bar dipping.
- Machine speed: up to 50 bars/minute (dry run).
- Machine output: up to 90,000 capsules/hour (depending on capsule size)

| Model No. | |
|------------------------|--|
| 7D2A-2-DECK-TES | Two-Deck Elongated Hard-Gelatin Machine |
| Specifications | |
| Dimensions (mm/inches) | 18,609.06mmL x 1651mmW x 2297.18mmH(mm) / 732.64" L x 65" W x 90.44" H (inches) |
| Weight (kg/lb) | Approx. (9797.59 kg / 21,600 lbs.) |

TECHNOPHAR

Technophar Global Headquarters

www.technophar.com

Updated: Dec. 2025
Copyright Technophar ©2025