

# Hard Gelatin Capsule Machine

7D2A 2-Deck Elongated Model



# Hard Gelatin Capsule Machine

## 7D2A 2-Deck Elongated Model

### Efficient Hard Gel Capsule Drying

The 7D2A 2-Deck TES elongated model has been designed to provide automated capsule filling capabilities for Hard Gel Capsules. Designed to provide high efficiency for production scale use meeting GMP standards for both pharmaceutical and nutraceutical industries.

#### The system consists of:

- Two deck machine:   Upper Deck: 30 tables  
                                  Lower Deck: 25 tables
- Control of all production and process parameters using high-tech control panel.
- Very easy and flexible operation using a Programmable Logic Controller (PLC) with touchscreen.

### Features:

- Two sets of automatics.
- Dipping 7-bars.
- 28 bars on table.
- Pinbar requirements: 3,208 pieces for 7-bar dipping.
- Machine speed: up to 50 bars/minute (dry run).
- Machine output: up to 90,000 capsules/hour (depending on capsule size)

Model No.	
7D2A-2-DECK-TES	Two-Deck Elongated Hard-Gelatin Machine
Specifications	
Dimensions (mm/inches)	18,609.06mmL x 1651mmW x 2297.18mmH(mm) / 732.64" L x 65" W x 90.44" H (inches)
Weight (kg/lb)	Approx. (9797.59 kg / 21,600lbs)

## **TECHNOPHAR**

Technophar Global Headquarters

[www.technophar.com](http://www.technophar.com)

Copyright Technophar ©2022 2022-01-18rev.